PULSE OXIMETER
Instruction manual

Reorder No：OXIM-L01-WL

Thank you for choosing WELL. Please read carefully the following instructions and keep them within reach.

- Batteries must be removed if the device is going to be stored for more than one month, or else batteries may leak.
- A flexible circuit connects the two parts of the device. Do not twist or pull on the connection.
- This pulse oximeter release audible warning sounds when pulse or oxygen values is nor within the standard.

1.4 Indicator for Use



The Fingertip Pulse Oximeter is a device intended for the spot-check of oxygen saturation of arterial hemoglobin (SpO²) and the pulse rate of adult. This device is not intended for continuous monitoring.

2. OVERVIEW

The pulse oxygen saturation is the percentage of HbO² in the total Hb in the blood, so-called the O² concentration in the blood. It is an important bio-parameter for the respiration. For the purpose of measuring the SpO² more easily and accurately, our company developed the Pulse Oximeter. At the same time, the device can measure the pulse rate simultaneously.

The Pulse Oximeter features in small volume, low power consumption, convenient operation and being portable. It is only necessary for test to put one of his fingers into a fingertip photoelectric sensor for diagnosis, and a display screen will directly show measured value of Hemoglobin Saturation.

Digital tube

6.2 Battery

- Step 1. Refer to Figure 3 and insert the two AAA size batteries properly in the right direction.
- Step 2. Replace the cover.
- ▲ Please take care when you insert the batteries for the improper insertion may damage the device.

INSTRUCTIONS TO USER

The Manual is written for the current Pulse Oximeter. In case of modifications and software upgrades, the information contained in this document is subject to change without notice.

The Manual describes in accordance with the Pulse Oximeter's features and requirements, main structure, functions, specifications, correct methods for transportation, installation, usage, operation, repair, maintenance and storage, etc.as well as the safety procedures to protect both the user and equipment. Refer to the respective chapters for details.

Please read the Manual very carefully before using this equipment.The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, personal injury and equipment damage due to user negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

WARNING

- The uncomfortable or painful feeling may appear if using the device ceaselessly. It is recommended that the sensor should not be applied to the same finger for over 2 hours.
- The light (the infrared is invisible) emitted from the

- device is harmful to the eyes, so the user and the maintenance man, cannot stare at the light.
- Test fingernail cannot be too long.
 - Please peruse the relative content about the clinical restrictions and caution.
 - This device is not intended for treatment.

1. SAFETY

1.1 Instructions for Safe Operations

- Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect.
 - It is recommended that the device should be inspected once a week at least. When there is obvious damage, stop using the monitor.
 - Necessary maintenance must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves.
 - The oximeter cannot be used together with devices not specified in User Manual. Only the accessory that appointed or recommendatory by manufacture can be used with this device.
- This product is calibrated before leaving factory.


1.2 Warnings

- Explosive hazard: Do NOT use the oximeter in environment with flammable gas such as some

3.1 Principle of Measurement

Principle of the Oximeter is as follows: An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive Hemoglobin (Hb) and Oxyhemoglobin (HbO²).

Operation principle of the instrument is: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning Recording Technology, so that two beams of different wavelength of lights can be focused onto human nail tip throughperspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on screen through treatment in electronic circuits and microprocessor.

8. REPAIRING AND MAINTENANCE

- Please change the batteries when the low-voltage displayed on the screen.
- Please clean the surface of the device before using.
- Wipe the device with medical alcohol first and then let it dry in air or clean it by dry clean fabric.
- Using the alcohol to disinfect the product after use, prevent from cross infection for next time use.
- Please take out the batteries if the oximeter is not in use for a long time.
- The best storage environment of the device is -40℃ to 60℃ ambient temperature and not higher than 95% relative humidity.

- ▲ High-pressure sterilization cannot be used on the device.Do not immerse the device in liquid.
- ▲ It is recommended that the device should be kept in a dry environment.
- ▲ Humidity may reduce the useful life of the device, or even damage it.

ignitable anesthetic agents.

- DO NOT use the oximeter while measured by MRI and CT.
- The person who is allergic to rubber cannot use this device.
- The disposal of scrap instrument and its accessories and packings (including battery, plastic bags, foams and paper boxes) should follow the local laws and regulations.
- Please check the packing before use to make sure the device and accessories are totally in accordance with the packing list, or else the device may have the possibility of working abnormally.

1.3 Attention

- Keep the oximeter away from dust, vibration, corrosive substances, explosive materials,high temperature and moisture.
- If the oximeter gets wet, please stop operating it.
- When it is carried from cold environment to warm or humid environment please do not use it immediately.
- DO NOT operate keys on front panel with sharp materials.
- High temperature or high pressure steam disinfection of the oximeter is not permitted. Refer to User Manual in the relative chapter for instructions of cleaning and disinfection.
- Do not immersed the oximeter in liquid. When is needs cleaning, please wipe its surface with medical alcohol by soft material Do not spray any liquid on the device directly.
- When cleaning the device with water, the temperature

3.2 Caution

1. The finger should be placed properly (see the attached illustration of this manual, Figure 5), or else it may cause inaccurate measurement.
 2. The SpO² sensor and photoelectric receiving tube should be arranged in a way with the subject article in a position there between.
 3. Make sure the optical path is free from any optical obstacles like rubberized fabric.
 4. Excessive ambient light may affect the measuring result. It includes fluorescent lamp, infrared heater, direct sunlight and etc.
 5. Strenuous action of the subject or extreme electrosurgical interference may also affect the accuracy.
- 3.3 Restrictions
1. As the measure is taken on the basis of arteriole pulse, substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO² waveform (PLETH)will decrease. In this case, the measurement will be more sensitive to interference.
 2. For those with a substantial amount of staining dilution drug (such as methylene blue), or carb monoxide hemoglobin (COHb), or methionine (Me + Hb) or Thiosalicylic hemoglobin, the SpO² determination by this monitor may be inaccurate.
 3. The drugs like dopamine, procaine, prilocaine,

- should be lower than 60℃.
- As to the fingers which are too thin or too cold, it would probably affect the normal measure of the SpO² test and pulse rate, please clip the thick finger such as thumb andmiddle finger deeply enough into the probe.
 - Do not use the device on infant.
 - The product is suitable for children above four years old and adults (Weight should be between 15 kg to 110 kg).
 - The device may not work in all tests. If you are unable to achieve stable readings, discontinue use.
 - The update period of data is less than 5 seconds, which is changeable according to different individual pulse rate.
 - The waveform is normalized. Please read the measured value when the waveform on screen is equably and steady-going, this measured value is optimal value, and the waveform at the moment is the standard one.
 - If some abnormal conditions appear on the screen during test process, pull out the finger and reinsert to restore normal use.
 - The device has normal useful life for three years since the first use.
 - The hanging rope attached the product is made from non-allergy material, if the personare sensitive to the hanging rope, stop using it. In addition, pay attention to the use of the hanging rope, do not wear it around the neck.
 - The instrument does not have low-voltage alarm function, it only shows the low-voltage. please change the battery when the battery energy is used out.

- lidocaine and butacaine may also be a major factor blamed for serious error of SpO² measure.
4. As the SpO² value serves as a reference value for judgement, someone with serious anemia may also report good SpO² measurement.

4. TECHNICAL SPECIFICATIONS

1. Type: Finger clip
2. Response Time:<5s
3. Battery: 2xAAA
4. Operating temperature: 5℃~40℃
5. Storage temperature: 10℃~50℃
6. Hemoglobin saturation:Upper limit: 100 / Bottom limit: 94
7. Pulse rate limit: Upper limit:130/ Bottom limit: 50
8. Hemoglobin saturation display: 35~100%
9. Pulse rate display:30~250BPM
10. Dimension:61.8X34.2X33.9MM
11. Net Weight:27.8g

5. ACCESSORIES

- One hanging rope.
- One User Manual.

6. INSTALLATION

6.1 Screen

10.THE MEANING OF SYMBOLS

Symbol	Description
	Type BF
	Refer to instruction manual/booklet
%SpO ²	The pulse oxygen saturation(%)
PRbpm	Pulse rate(bpm)
	The battery voltage indication is deficient (change the battery in time avoiding the inexact measure)
	1, No finger inserted 2, An indicator of signal inadequacy
+	Battery positive e lectrode
-	Battery cathode
	Enterprise related information
	Authorized repres entative in the European community
	The symbol indicates that the device complieswith the European Council Directive 93/42/EEC concerning medical devices
	WEEE(2002/96/EC)
IP22	International Protection

11. FUNCTION SPECIFICATION

Display Information	Display Mode
The Pulse Oxygen Saturation (SpO ²)	Digital display
Pulse Rate (PR)	Digital display
Pulse Intensity (bar-graph)	Bar-graph display
SpO ² Parameter Specification	
Measuring range	1% ~100%, (resolution is 1%)
Accuracy	70% ~100%, (resolution is 2%)
Optical sensor	Red light (wavelength is 660nm) Infrared (wavelength is 880nm)
Pulse Parameter Specification	
Measuring range	25 bpm 250 bpm (resolution is 1 bpm)
Accuracy	± 2 bpm or ± 2%
Pulse Intensity	
Range	Continuous bar-graph display, the higher display indicates the strong pulse
Battery Requirement	
1.5 V (AAA size) alkaline batteries x 2 or rechargeable battery	

Battery Useful Life	
Two batteries can work continually for 20 hours	
Dimensions and Weight OXIM-L01-WL	
Dimensions	61,8X34,2X33,9MM
Weight	About 27.8g

	Name: MedPath G mbH Address: Mies-van-der-Rohe-Strasse 8,80807 Munich, Germany Telephone: +49(0)89 189174474 E-Mail: info@medpathpro
--	---

Waste electrical and electronic equipment are a special waste category, collection, storage, transport, treatment and recycling are important because they can avoid environmental pollution and are harmful to health.

Submitting waste electrical and electronic equipment to special collection centers makes the waste to be recycled properly and protecting the environment.

Do not forget! Each electric appliance that arrives at the landfill, the field, pollutes the environment!

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